

NEB #	DNA Assembly, Cloning and Mutagenesis Kits	Size	Stockroom #
E2621S	NEBuilder® HiFi DNA Assembly Master Mix	10 reactions	E2621S
E0554S	Q5® Site-Directed Mutagenesis Kit	10 reactions	E0554S
NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M0525S	Quick CIP	1,000 units	M0525S
M2200S	Quick Ligation™ Kit	30 reactions	M2200S
M0202S	T4 DNA Ligase	20,000 units	M0202S
M0202L	T4 DNA Ligase	100,000 units	M0202L
M0201S	T4 Polynucleotide Kinase	500 units	M0201S
M0315S	Terminal Transferase	500 units	M0315S
NEB #	Markers & Ladders	Size	Stockroom #
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	N3232S
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	N3231S
N0467S	Quick-Load® 100 bp DNA Ladder	125 gel lanes (1.25 ml)	N0467S
NEB #	Nucleic Acid Purification	Size	Stockroom #
T2010S	Monarch® Total RNA Miniprep Kit	50 preps	T2010S
T1120S	Monarch® Spin DNA Gel Extraction Kit	50 preps	T1120S
T1130S	Monarch® Spin PCR & DNA Cleanup Kit (5 µg)	50 preps	T1130S
T1130L	Monarch® Spin PCR & DNA Cleanup Kit (5 µg)	250 preps	T1130L
T1110S	Monarch® Spin Plasmid Miniprep Kit	50 preps	T1110S
T1110L	Monarch® Spin Plasmid Miniprep Kit	250 preps	T1110L
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
M0538S	Bst 2.0 WarmStart® DNA Polymerase	1,600 units	M0538S
N0447S	Deoxynucleotide (dNTP) Solution Mix	8 µmol	N0447S
M0480S	OneTaq® DNA Polymerase	200 units	M0480S
M0481S	OneTaq® Hot Start DNA Polymerase	200 units	M0481S
M0530S	Phusion™ High-Fidelity DNA Polymerase	100 units	M0530S
M0530L	Phusion™ High-Fidelity DNA Polymerase	500 units	M0530L
M0491S	Q5® High-Fidelity DNA Polymerase	100 units	M0491S
M0491L	Q5® High-Fidelity DNA Polymerase	500 units	M0491L
M0273S	Taq DNA Polymerase with Standard Taq Buffer	400 units	M0273S
M0267S	Taq DNA Polymerase with ThermoPol® Buffer	400 units	M0267S
M0267L	Taq DNA Polymerase with ThermoPol® Buffer	2,000 units	M0267L

NEB #	Protein Expression & Purification Technologies	Size	Stockroom #
S6651S	Chitin Resin	20 ml	S6651S
NEB #	Restriction Endonucleases	Size	Stockroom #
R3552S	Agel-HF®	300 units	R3552S
R3136S	BamHI-HF®	10,000 units	R3136S
R0176S	Dpnl	1,000 units	R0176S
R3101S	EcoRI-HF®	10,000 units	R3101S
R3195S	EcoRV-HF®	4,000 units	R3195S
R3104S	HindIII-HF®	10,000 units	R3104S
R3142S	KpnI-HF®	4,000 units	R3142S
R3193S	NcoI-HF®	1,000 units	R3193S
R0111S	NdeI	4,000 units	R0111S
R3131S	NheI-HF®	1,000 units	R3131S
R3189L	NotI-HF®	2,500 units	R3189L
R3189S	NotI-HF®	500 units	R3189S
R0547S	PacI	250 units	R0547S
R3140S	PstI-HF®	10,000 units	R3140S
R3156S	SacI-HF®	2,000 units	R3156S
R3138S	Sall-HF®	2,000 units	R3138S
R3122S	ScaI-HF®	1,000 units	R3122S
R0141S	SmaI	2,000 units	R0141S
R3182S	SphI-HF®	500 units	R3182S
R0145S	XbaI	3,000 units	R0145S
R0146S	XhoI	5,000 units	R0146S
R0180S	XmaI	500 units	R0180S

NEB Representative: Shanna Sawatzki
shanna@neb.com
(206) 889-0401

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com