

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M2200S	Quick Ligation™ Kit	30 reactions	M2200S
M0371S	Shrimp Alkaline Phosphatase (rSAP)	500 units	M0371S
M0202S	T4 DNA Ligase	20,000 units	M0202S
M0201S	T4 Polynucleotide Kinase	500 units	M0201S
NEB #	Glycobiology	Size	Stockroom #
P0704S	PNGase F	15,000 units	P0704S
NEB #	Markers & Ladders	Size	Stockroom #
N3232L	1 kb DNA Ladder	1,000 gel lanes (1 ml)	N3232L
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	N3232S
N3200S	1 kb Plus DNA Ladder	200 gel lanes (0.1 ml)	N3200S
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	N3231S
N0467S	Quick-Load® 100 bp DNA Ladder	125 gel lanes (1.25 ml)	N0467S
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
M0258S	Deep Vent® DNA Polymerase	200 units	M0258S
N0447S	Deoxynucleotide (dNTP) Solution Mix	8 µmol	N0447S
M0210S	DNA Polymerase I, Large (Klenow) Fragment	200 units	M0210S
M0484S	OneTaq® Hot Start 2X Master Mix with Standard Buffer	100 reactions (50 µl)	M0484S
M0530S	Phusion™ High-Fidelity DNA Polymerase	100 units	M0530S
M0368S	ProtoScript® II Reverse Transcriptase	4,000 units	M0368S
M0491S	Q5® High-Fidelity DNA Polymerase	100 units	M0491S
M0320S	Taq DNA Polymerase with Standard Taq (Mg-free) Buffer	400 units	M0320S
NEB #	Restriction Endonucleases	Size	Stockroom #
R0136S	BamHI	10,000 units	R0136S
R0176S	DpnI	1,000 units	R0176S
R0101S	EcoRI	10,000 units	R0101S
R0195S	EcoRV	4,000 units	R0195S
R0193S	NcoI	1,000 units	R0193S
R0189S	NotI	500 units	R0189S
R0140S	PstI	10,000 units	R0140S
R0138S	SalI	2,000 units	R0138S
R0141S	SmaI	2,000 units	R0141S
R0182S	SphI	500 units	R0182S

R0146S	Xhol	5,000 units	R0146S
--------	------	-------------	--------

NEB Representative: Jens Erichsen
jerichsen@neb.com
(919) 306-5444

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com