

B7024S	Gel Loading Dye, Purple (6X)	4 ml	1
B9200S	Recombinant Albumin, Molecular Biology Grade	12 mg	1
E2621S	NEBuilder® HiFi DNA Assembly Master Mix	10 reactions	1
M0289S	Antarctic Phosphatase	1,000 units	2
M0525S	Quick CIP	1,000 units	1
M2200S	Quick Ligation™ Kit	30 reactions	1
M0202S	T4 DNA Ligase	20,000 units	2
M0202L	T4 DNA Ligase	100,000 units	1
N0552S	Quick-Load® Purple 1 kb DNA Ladder	125 gel lanes (1.25 ml)	1
N0551L	Quick-Load® Purple 100 bp DNA Ladder	375 gel lanes (3.75 ml)	1
T3050L	Monarch® HMW DNA Extraction Kit for Cells & Blood	50 preps	1
T3060L	Monarch® HMW DNA Extraction Kit for Tissue	50 preps	1
T1030L	Monarch® PCR & DNA Cleanup Kit (5 µg)	250 preps	1
T3010S	Monarch® Spin gDNA Extraction Kit	50 preps	1
N0447S	Deoxynucleotide (dNTP) Solution Mix	8 µmol	1
M0210L	DNA Polymerase I, Large (Klenow) Fragment	1,000 units	1
M0323L	LongAmp® Taq DNA Polymerase	2,500 units	2
M0323S	LongAmp® Taq DNA Polymerase	500 units	1
M0541S	NEBNext® High-Fidelity 2X PCR Master Mix	50 reactions	2
M0480S	OneTaq® DNA Polymerase	200 units	1
M0530L	Phusion® High-Fidelity DNA Polymerase	500 units	1
M0530S	Phusion® High-Fidelity DNA Polymerase	100 units	1
M0491S	Q5® High-Fidelity DNA Polymerase	100 units	2
M0491L	Q5® High-Fidelity DNA Polymerase	500 units	2
M0494S	Q5® Hot Start High-Fidelity 2X Master Mix	100 reactions	1
R3552S	AgeI-HF®	300 units	1
R3136S	BamHI-HF®	10,000 units	2
R3101S	EcoRI-HF®	10,000 units	1
R3104S	HindIII-HF®	10,000 units	2
R3193S	NcoI-HF®	1,000 units	1
R3189S	NotI-HF®	500 units	1