

NEB #	Buffers	Size	Quantity Available
B7024S	Gel Loading Dye, Purple (6X)	4 ml	2
NEB #	DNA Assembly, Cloning and Mutagenesis Kits	Size	Quantity Available
E2611S	Gibson Assembly® Master Mx	10 reactions	2
E2621S	NEBuilder® HiFi DNA Assembly Master Mx	10 reactions	2
E0552S	Q5® Site-Directed Mutagenesis Kit (Without Competent Cells)	10 reactions	2
NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Quantity Available
M0303S	DNase I (RNase-free)	1,000 units	1
M0525S	Quick CIP	1,000 units	2
M2200L	Quick Ligation™ Kit	150 reactions	1
M2200S	Quick Ligation™ Kit	30 reactions	2
M0371S	Shrimp Alkaline Phosphatase (rSAP)	500 units	1
M0202T	T4 DNA Ligase	20,000 units	2
M0202L	T4 DNA Ligase	100,000 units	2
M0202S	T4 DNA Ligase	20,000 units	3
B0202S	T4 DNA Ligase Reaction Buffer	5 ml	2
M0201S	T4 Polynucleotide Kinase	500 units	3
NEB #	Markers & Ladders	Size	Quantity Available
N3232L	1 kb DNA Ladder	1,000 gel lanes (1 ml)	1
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	out of stock
N3200L	1 kb Plus DNA Ladder	1,000 gel lanes (0.5 ml)	1
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	out of stock
N3231L	100 bp DNA Ladder	500 gel lanes (0.5 ml)	out of stock
N0468S	Quick-Load® 1 kb DNA Ladder	125 gel lanes (1.25 ml)	1
N0467S	Quick-Load® 100 bp DNA Ladder	125 gel lanes (1.25 ml)	2
N0550S	Quick-Load® Purple 1 kb Plus DNA Ladder	250 gel lanes (1.25 ml)	2
N0551S	Quick-Load® Purple 100 bp DNA Ladder	125 gel lanes (1.25 ml)	2
NEB #	Nucleic Acid Purification	Size	Quantity Available
T2010S	Monarch® Total RNA Miniprep Kit	50 preps	1
E1050S	Exo-CIP™ Rapid PCR Cleanup Kit	100 reactions	2
T2040L	Monarch® RNA Cleanup Kit (50 µg)	100 preps	1
T2040S	Monarch® RNA Cleanup Kit (50 µg)	10 preps	1
T1020L	Monarch® DNA Gel Extraction Kit	250 preps	1
T1030L	Monarch® PCR & DNA Cleanup Kit (5 µg)	250 preps	1
T1030S	Monarch® PCR & DNA Cleanup Kit (5 µg)	50 preps	1
T1010S	Monarch® Plasmid Miniprep Kit	50 preps	1
T1010L	Monarch® Plasmid Miniprep Kit	250 preps	1
T1120S	Monarch® Spin DNA Gel Extraction Kit	50 preps	1
T3010S	Monarch® Spin gDNA Extraction Kit	50 preps	1
T3010L	Monarch® Spin gDNA Extraction Kit	150 preps	1
T1130S	Monarch® Spin PCR & DNA Cleanup Kit (5 µg)	50 preps	1

NEB #	PCR, Polymerases & Amplification Technologies	Size	Quantity Available
N0447S	Deoxynucleotide (dNTP) Solution Mix	8 µmol	1
N0446S	Deoxynucleotide (dNTP) Solution Set	25 µmol	2
E3005S	Luna [®] Universal One-Step RT-qPCR Kit	200 reactions	2
M8003S	Luna [®] Universal qPCR Master Mix	200 reactions	2
M0481S	OneTaq [®] Hot Start DNA Polymerase	200 units	2
M0488S	OneTaq [®] Hot Start Quick-Load [®] 2X Master Mix with Standard Buffer	100 reactions (50 µl)	2
M0486S	OneTaq [®] Quick-Load [®] 2X Master Mix with Standard Buffer	100 reactions (50 µl)	2
M0530S	Phusion [®] High-Fidelity DNA Polymerase	100 units	2
E6560S	ProtoScript [®] II First Strand cDNA Synthesis Kit	30 reactions	1
M0368S	ProtoScript [®] II Reverse Transcriptase	4,000 units	1
M0491L	Q5 [®] High-Fidelity DNA Polymerase	500 units	1
M0491S	Q5 [®] High-Fidelity DNA Polymerase	100 units	1
M0494S	Q5 [®] Hot Start High-Fidelity 2X Master Mix	100 reactions	1
M0493S	Q5 [®] Hot Start High-Fidelity DNA Polymerase	100 units	2
M0297S	RNase H	250 units	1
M0314S	RNase Inhibitor, Murine	3,000 units	2
M0203S	T4 DNA Polymerase	150 units	1
NEB #	Protein Tools	Size	Quantity Available
P8107S	Proteinase K, Molecular Biology Grade	2 ml	2
NEB #	Restriction Endonucleases	Size	Quantity Available
R3552S	AgeI-HF [®]	300 units	3
R0558S	Ascl	500 units	2
R3136S	BamHI-HF [®]	10,000 units	2
R0144S	BglII	2,000 units	2
R0176S	DpnI	1,000 units	2
R3101S	EcoRI-HF [®]	10,000 units	2
R3195S	EcoRV-HF [®]	4,000 units	2
R3104S	HindIII-HF [®]	10,000 units	2
R3142S	KpnI-HF [®]	4,000 units	2
R3193S	NcoI-HF [®]	1,000 units	2
R3131S	NheI-HF [®]	1,000 units	2
R3189S	NotI-HF [®]	500 units	2
R3189L	NotI-HF [®]	2,500 units	2
R0547S	Pacl	250 units	2
R3156S	SacI-HF [®]	2,000 units	2
R3138S	Sall-HF [®]	2,000 units	2
R0141S	SmaI	2,000 units	1
R0145S	XbaI	3,000 units	2
R0146S	XhoI	5,000 units	2

NEB Representative: Lynn McBee
 mcbee@neb.com
 (214) 707-7065

NEBnow Support: 1-800-632-5227
 Email: freezers@neb.com